



SIEVE SHAKER OCTAGON 200CL AND SIEVEWARE EVALUATION SOFTWARE

Endecotts sieve shaker Octagon 200CL competes with the most advanced sieve shakers in the world. Several unique features have been developed specifically for this machine, including the 'Closed Loop' amplitude control for ultimate reproducibility.

It is also fully compatible with the SieveWare - evaluation and control software via RS232 Port.

Endecotts evaluation software makes sieve analysis simple and more comfortable than manual evaluation.

The SieveWare enables the automatic sieve evaluation with the Octagon 200CL.

Once set up, the evaluation and weighing process becomes extremely easy and efficient. It eliminates errors and provides trouble-free analysis of your sieving data.





SIEVE SHAKER OCTAGON 200CL

For precise, reproducible and error-free sieving processes



Easy-to-use
sieve clamping

Accepts up to
8 full height 200 mm / 8"
diameter sieves


Suitable for dry and
wet sieving

3D sieving action
allows for high
separation efficiency

'Closed Loop' total
amplitude control
ensures reproducible
sieving

Digital controls for
easy and reliable
operation

RS232 Port



Full compatibility
with SieveWare
evaluation software

Endecotts Ltd
Parsons Lane, Hope
Hope Valley, S33 6RB
United Kingdom

Telephone: +44 (0)20 8542 8121
E-mail: sales@endecotts.com
Web: www.endecotts.com

Endecotts US Customer Service
Toll free number: 855 897 7284
E-mail: customer.service@endecotts-us.com

WWW.ENDECOTTS.COM